


Proper Fit Every Time

It is crucial to measure stoma length using the AMT Balloon Stoma Measuring Device for proper fit, especially when a patient switches to AMT from another brand.

Balloon Stoma Measuring Device	Order #	Packaged
Non-Guidewire Compatible	1704	10/Box
Guidewire Compatible	1804	10/Box



F = Tube Diameter
cm = Stoma Length

Example: 14F x 2.3 cm

(F) French Size
(cm) Centimeter

MiniONE® Button

Low Profile Balloon Gastrostomy Tube

MiniONE® Balloon Button Sizes

12 F		14 F		16 F		18 F		20 F		24 F	
(cm) Length	Order #	(cm) Length	Order #	(cm) Length	Order #	(cm) Length	Order #	(cm) Length	Order #	(cm) Length	Order #
0.8	M1-5-1208	0.8	M1-5-1408								
1.0	M1-5-1210	1.0	M1-5-1410	1.0	M1-5-1610	1.0	M1-5-1810	1.0	M1-5-2010		
1.2	M1-5-1212	1.2	M1-5-1412	1.2	M1-5-1612	1.2	M1-5-1812	1.2	M1-5-2012		
1.5	M1-5-1215	1.5	M1-5-1415	1.5	M1-5-1615	1.5	M1-5-1815	1.5	M1-5-2015	1.5	M1-5-2415
1.7	M1-5-1217	1.7	M1-5-1417	1.7	M1-5-1617	1.7	M1-5-1817	1.7	M1-5-2017	1.7	M1-5-2417
2.0	M1-5-1220	2.0	M1-5-1420	2.0	M1-5-1620	2.0	M1-5-1820	2.0	M1-5-2020	2.0	M1-5-2420
2.3	M1-5-1223	2.3	M1-5-1423	2.3	M1-5-1623	2.3	M1-5-1823	2.3	M1-5-2023	2.3	M1-5-2423
2.5	M1-5-1225	2.5	M1-5-1425	2.5	M1-5-1625	2.5	M1-5-1825	2.5	M1-5-2025	2.5	M1-5-2425
2.7	M1-5-1227	2.7	M1-5-1427	2.7	M1-5-1627	2.7	M1-5-1827	2.7	M1-5-2027	2.7	M1-5-2427
3.0	M1-5-1230	3.0	M1-5-1430	3.0	M1-5-1630	3.0	M1-5-1830	3.0	M1-5-2030	3.0	M1-5-2430
3.5	M1-5-1235	3.5	M1-5-1435	3.5	M1-5-1635	3.5	M1-5-1835	3.5	M1-5-2035	3.5	M1-5-2435
4.0	M1-5-1240	4.0	M1-5-1440	4.0	M1-5-1640	4.0	M1-5-1840	4.0	M1-5-2040	4.0	M1-5-2440
		4.4	M1-5-1444	4.4	M1-5-1644	4.4	M1-5-1844	4.4	M1-5-2044	4.4	M1-5-2444
		5.0	M1-5-1450			5.0	M1-5-1850				
		5.5	M1-5-1455								
		6.0	M1-5-1460			6.0	M1-5-1860				
		6.5	M1-5-1465								

All MiniONE kits above are packaged with one 8-1255G Y-Port feed set and one 8-1211 Bolus feed set.

Compatible Feed Sets

		Connector Options							
		Right Angle (GG)		Right Angle		Straight			
		12"	24"	2"	12"	24"	2"	12"	24"
Adapter Options	Y-Port (Purple)	8-1255G*		8-1255	8-2455				
	Bolus	8-1257G**		8-1221	8-2421	8-0211	8-1211*	8-2411	
	Y-Port (Clear)			8-0222	8-1222	8-2422	8-0212	8-1212	

GG=Glow Green

*Included in MiniONE® Kits
**Y-Port with Enhanced Med Port Plug
All come in a package of 10/box.

All MiniONE® feed sets are compatible with Mic-Key® devices.

Note: Feeding sets are supplied assembled. Choose the best combination of length, connector and adapter type for your use. Contact your insurance provider regarding the number of feeding/ extension sets allowed per month.

AMTClamp™

Order # 4-3000



For a complete list of products visit us at www.AppliedMedical.net

800 869 7382

CS@AppliedMedical.net



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net

An ISO 13485 Registered Company CE

© 2016 Applied Medical Technology, Inc.



US Patent Numbers 6,019,746, other patents pending. Mic-Key® is a registered trademark or trademark of Halyard Worldwide., or its affiliates C-3922-G



Reduce Leaks & Improve Comfort

Improve Patient Outcomes

Easy Nighttime use Glows in the Dark



Innovating. Educating. Changing Lives.



What Makes A Better Balloon Button?

Reduce Leaks


The AMT MiniONE® is designed to maximize seal at the gastric surface for a better fit.

- **Apple Shaped Balloon** grows horizontally when filled, covering more surface area at the gastrostomy site.
- **Tapered Transition** from the shaft to the base of the external bolster is designed for a more secure fit.
- **Durable** anti-reflux valve offers maximum protection.
- **Adjustable** balloon allows for fill variation to ensure individualized proper fit.



The AMT MiniONE® is designed to best fulfill a patient's nutritional needs with as little lifestyle disruption as possible.

37% Lower Profile External Bolster



40% Lower Profile Within The Stomach

Balloon Fill Volumes

French Size	Minimum	Recommended	Maximum
12F	2 ml	2.5 ml	3 ml
14F	3 ml	4 ml	5 ml
16F	4 ml	6 ml	8 ml
18F	6 ml	8 ml	10 ml
20F	7 ml	10 ml	15 ml
24F	7 ml	10 ml	15 ml

Consult your healthcare provider before adjusting fill volumes.

Recommended amount of water to fill the balloon is printed on fill valve



Improve Comfort

Maximization of comfort is key when it comes to patient, caregiver and clinician satisfaction with an enteral feeding device.

- **37% Lower* Profile** external bolster lays discreetly beneath clothing.
- **40% Lower* Profile** internal retention balloon occupies less gastric space.
- **Soft** flexible high grade silicone is friendly to skin surface and easy to clean around.
- **Unobtrusive Fill Port** angled away from the body for reduced irritation to skin.
- **Glow Green** feeding port and feed set connector glows in the dark to assist caregivers with administration of nighttime feeds avoiding the need to wake the patient.




Clearly Labeled



Glow Green for Nighttime Use



Softer



For Maximum Balloon Life

It is crucial to measure stoma length using the AMT Balloon Stoma Measuring Device for proper fit.



Guidewire Compatible
Order #: 1804



Non-Guidewire Compatible
Order #: 1704

AMT Balloon Stoma Measuring Device

*Lower than the leading competitor