


Dimensional Comparison

Dimensional Analysis: MiniONE® vs. Competitor



The Importance of Measuring

The “length” of the MiniONE refers to the distance from the bottom of the external bolster to the top of the balloon. This measurement may be slightly different with other brands.

Stoma length measurement is always important, but it is especially crucial to use the **AMT Balloon Stoma Measuring Device** when a patient switches to AMT from another brand. Due to differences in design, it cannot be assumed a patient will use the same size in both devices.

Contact Us:
800 869 7382
CS@AppliedMedical.net



US Patent Numbers 6,019,746 and other Patents Pending.

Applied Medical Technology, Inc.

www.AppliedMedical.net

Phone: 800 869 7382


Fax: 440 717 4200

CS@AppliedMedical.net

C-4495-A
Rev 1/16

Balloon Volume Comparison

Gastric Seal Surface Analysis: MiniONE® vs. Competitor



Greater Comfort

The MiniONE® Balloon Button is designed to best fulfill nutritional needs with as little lifestyle disruption as possible.



“Apple”-shaped balloon grows horizontally when filled instead of vertically.