



Patient Name _____

DO NOT SUBSTITUTE

Refill Allowed



ENFit® Kit



Legacy Kit

**MiniONE® Button Low Profile Balloon Gastrostomy Tube**

12 French			14 French			16 French			18 French			20 French			24 French		
Length	Order #		Length	Order #		Length	Order #		Length	Order #		Length	Order #		Length	Order #	
0.8	M1-5-1208		0.8	M1-5-1408													
1.0	M1-5-1210		1.0	M1-5-1410		1.0	M1-5-1610		1.0	M1-5-1810		1.0	M1-5-2010				
1.2	M1-5-1212		1.2	M1-5-1412		1.2	M1-5-1612		1.2	M1-5-1812		1.2	M1-5-2012				
1.5	M1-5-1215		1.5	M1-5-1415		1.5	M1-5-1615		1.5	M1-5-1815		1.5	M1-5-2015		1.5	M1-5-2415	
1.7	M1-5-1217		1.7	M1-5-1417		1.7	M1-5-1617		1.7	M1-5-1817		1.7	M1-5-2017		1.7	M1-5-2417	
2.0	M1-5-1220		2.0	M1-5-1420		2.0	M1-5-1620		2.0	M1-5-1820		2.0	M1-5-2020		2.0	M1-5-2420	
2.3	M1-5-1223		2.3	M1-5-1423		2.3	M1-5-1623		2.3	M1-5-1823		2.3	M1-5-2023		2.3	M1-5-2423	
2.5	M1-5-1225		2.5	M1-5-1425		2.5	M1-5-1625		2.5	M1-5-1825		2.5	M1-5-2025		2.5	M1-5-2425	
2.7	M1-5-1227		2.7	M1-5-1427		2.7	M1-5-1627		2.7	M1-5-1827		2.7	M1-5-2027		2.7	M1-5-2427	
3.0	M1-5-1230		3.0	M1-5-1430		3.0	M1-5-1630		3.0	M1-5-1830		3.0	M1-5-2030		3.0	M1-5-2430	
3.5	M1-5-1235		3.5	M1-5-1435		3.5	M1-5-1635		3.5	M1-5-1835		3.5	M1-5-2035		3.5	M1-5-2435	
4.0	M1-5-1240		4.0	M1-5-1440		4.0	M1-5-1640		4.0	M1-5-1840		4.0	M1-5-2040		4.0	M1-5-2440	
			4.4	M1-5-1444		4.4	M1-5-1644		4.4	M1-5-1844		4.4	M1-5-2044		4.4	M1-5-2444	
			5.0	M1-5-1450					5.0	M1-5-1850							
			5.5	M1-5-1455													
			6.0	M1-5-1460					6.0	M1-5-1860							
			6.5	M1-5-1465													

Balloon Stoma Measuring Device

Description	Order #
Non-Guidewire Compatible	1704
Guidewire Compatible	1804

See Reverse For Feed Sets.

MiniONE® ENFit® kits are available! – Order numbers will be the same as above with “-I” as a suffix for ENFit®, LEGACY MiniONE® KITS are packaged with one 8-1255 Y-Port feed set and one 8-1211 Bolus feed set. ENFit® MiniONE® KITS are packaged with one 8-1255-ISOSAF and one 8-1211-ISOSAF feed set, along with two ENFit® syringes.

MD Signature _____

Date _____

For instructional videos and more information: www.AppliedMedical.net

©2017 Applied Medical Technology, Inc. C-4339-F

Patents: <http://www.AppliedMedical.net/Patents/> – Trademarks: <http://www.AppliedMedical.net/Trademarks/>

MR Conditional

MiniONE® Feed Set Options

Order #	Length	Connector GG	Adapter
8-1255	12"	Right Angle	Purple Y-Port
8-2455	24"	Right Angle	Purple Y-Port
8-1257G*	12"	Right Angle	Purple Y-Port
8-1221	12"	Right Angle	Bolus
8-2421	24"	Right Angle	Bolus
8-0211	2"	Straight	Bolus
8-1211	12"	Straight	Bolus
8-2411	24"	Straight	Bolus
8-0222	2"	Right Angle	Clear Y-Port
8-1222	12"	Right Angle	Clear Y-Port
8-2422	24"	Right Angle	Clear Y-Port
8-0212	2"	Straight	Clear Y-Port
8-1212	12"	Straight	Clear Y-Port
8-1255-H	12"	Right Angle	Hybrid Y-Port
8-2455-H	24"	Right Angle	Hybrid Y-Port
8-1255-ISOSAF	12"	Right Angle	ENFit® Y-Port
8-2455-ISOSAF	24"	Right Angle	ENFit® Y-Port
8-1221-ISOSAF	12"	Right Angle	ENFit® Single Port
8-2421-ISOSAF	24"	Right Angle	ENFit® Single Port
8-0211-ISOSAF	02"	Straight	ENFit® Single Port
8-1211-ISOSAF	12"	Straight	ENFit® Single Port
8-2411-ISOSAF	24"	Straight	ENFit® Single Port

GG=Glow Green

*Included in MiniONE® Legacy Kits

*Included in MiniONE® ENFit® Kits

*Y-Port with Enhanced Med Port Plug. All come in a package of 10/box.

All MiniONE® feed sets are compatible with Mic-Key® devices.

Note: Feed sets are supplied assembled. Choose the best combination of length, connector and adapter type for your use. Contact your insurance provider regarding the number of Feed/extension sets allowed per month.




AMT CLAMP™

Order #	Description
<input type="checkbox"/>	For Use With Enteral Feed Tubes, Adapters, Y-Ports and G-Tubes
4-3000	

Connector Options

COMPATIBLE FEED SETS		Right Angle (GG)			Straight (GG)		
		2"	12"	24"	2"	12"	24"
Adapter Options	Purple Y-Port		8-1255 8-1257G*	8-2455			
	Bolus		8-1221	8-2421	8-0211	8-1211	8-2411
	Clear Y-Port	8-0222	8-1222	8-2422	8-0212	8-1212	
	Hybrid Y-Port		8-1255-H	8-2455-H			
	ENFit® Y-Port		8-1255-ISOSAF	8-2455-ISOSAF			
	ENFit® Single Port		8-1221-ISOSAF	8-2421-ISOSAF	8-0211-ISOSAF	8-1211-ISOSAF	8-2411-ISOSAF

Inflated Balloon



Deflated Balloon



F = Tube Diameter
cm = Gastric Length

For a complete list of products and accessories visit

www.AppliedMedical.net